

Antigen Preparation Prior to Injection

- 1.) Make 1ml of final emulsion. This is enough for two boosts. Freeze at -20°C after the first boost, and reemulsify before injecting again.
- 2.) Keep emulsifying needles at 4°C well labeled, and reuse the same one for the same antigen so we don't have to clean each time.
- 3.) Make antigen/adjuvant etc. to mix in 3ml syringe (luer-lok). Use pipetteman to add antigen into syringe "nose" after adding Freund's adjuvant by sucking up into syringe through a 20g or larger bore needle.
- 4.) When using Complete Freund's adjuvant (Gibco), warm first and make sure BCG are suspended. When using Incomplete Freund's, warm it so it sucks well through a needle.
- 5.) Inject approximately $100\mu\text{g}$ KLH-linked peptide for each injection. $200\mu\text{g}$ aliquots of KLH-linked peptide are stored at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ in $500\mu\text{l}$ size "O-ring" microcentrifuge tubes. $200\mu\text{g}$ aliquots have $\sim 200-400\mu\text{l}$ liquid.
- 6.) Mix - Antigen ($\sim 300-400\mu\text{l}$)
 PBS up to $500\mu\text{l}$ total
 Freunds 500ml or a little more if need to emulsify well

 Total $\sim 1\text{ml}$.
- 7.) In the past, we have injected $\sim 250\mu\text{l}$ / site, 2 sites per time. New Emory guidelines say $100\mu\text{l}$ of emulsion/site for rabbits at Emory.